


CHOICE STRETCH

HORIZONTAL GEL SYSTEMS



The Clarit-E Choice Stretch tanks are ideal for researchers wanting to perform higher resolution separation of more samples over a longer distance. The extended run lengths allow restriction fragments or other close molecular weight sample bands to be easily separated and identified. Speed loading is accomplished using 10, 14, 16, 18, 28 or 30 sample multichannel pipette compatible combs.

KEY FEATURES:

The Clarit-E Choice Stretch is deal for higher resolution separations such as separating samples with very similar molecular weights.

- α Available with 15 x 20cm and 15 x 25cm gel trays
- α Run up to 350 samples
- α Multichannel pipette compatible combs for speed loading



Casting Dams
Cast gels separately from the gel tank and save time while the tank is in use.



Buffer Saver Blocks
Save costs by physically reducing the volume of buffer required in the gel chamber.

Tanks:

| Product Code | Description |
|--------------|------------------|
| EL1430 | Gel Tank 15x20cm |
| EL1440 | Gel Tank 15x25cm |

Trays:

| Product Code | Description |
|--------------|-----------------|
| EL1432 | UV Tray 15x20cm |
| EL1442 | UV Tray 15x25cm |

Accessories:

| Product Code | Description |
|--------------|------------------------|
| EL1404 | Casting Dams |
| EL1406 | Electrode, positive |
| EL1408 | Electrode, negative |
| EL1410 | Loading Guides |
| EL1428 | Viewing Platform |
| EL1414 | Cool Pack and Platform |
| EL1416 | Buffer Saver Blocks |
| EL1418 | Flexi Caster |
| EL1026 | Gel Levelling Table |
| EL1522 | Lid |

Clarit-E®

Compatible Power Supplies and Adapters:

| Product Code | Description |
|--------------|--|
| EL2707 | Mini Power Supply 500V, 400mA, 120W |
| EL2705 | Mini-300 Power Supply 300V, 400mA, 60W |
| EL2702 | Power Pack Adapters 4mm to 2mm |
| EL2752 | Maxi Power Pack 300V, 3A, 300W |
| EL2732 | Midi Power Pack 300V, 700mA, 150W |
| EL2714 | Replacement Cables |

Compatible Combs:

| Product Code | Description |
|------------------|-------------------------------|
| EL1400-1-0.75 | Comb Prep 1, Marker 1, 0.75mm |
| EL1400-1-1 | Comb Prep 1, Marker 1, 1mm |
| EL1400-1-1.5 | Comb Prep 1, Marker 1, 1.5mm |
| EL1400-1-2 | Comb Prep 1, Marker 1, 2mm |
| EL1400-2-0.75 | Comb Prep 2, Marker 2, 0.75mm |
| EL1400-2-1 | Comb Prep 2, 1.0mm |
| EL1400-2-1.5 | Comb Prep 2, Marker 2, 1.5mm |
| EL1400-2-2 | Comb Prep 2, Marker 2, 2mm |
| EL1400-4-0.75 | Comb Prep 4, marker 2, 0.75mm |
| EL1400-4-1 | Comb Prep 4, marker 2, 1mm |
| EL1400-4-1-5 | Comb Prep 4, marker 2, 1.5mm |
| EL1400-4-2 | Comb Prep 4, marker 2, 2mm |
| EL1400-10-0.75 | Comb 10 Sample, 0.75mm |
| EL1400-10-1 | Comb 10 well, 1mm thick |
| EL1400-10-1.5 | Comb 10 well, 1.5mm |
| EL1400-10-2 | Comb 10 well, 2mm |
| EL1400-10MC-0.75 | Comb 10 Sample MC, 0.75 mm |
| EL1400-10MC-1 | Comb 10 well MC, 1mm |
| EL1400-10MC-1.5 | Comb 10 well MC, 1.5mm |
| EL1400-10MC-2 | Comb 12 well MC, 2mm |
| EL1400-12-0.75 | Comb 10 well MC, 0.75 mm |
| EL1400-12-1 | Comb 12 Sample, 1mm |
| EL1400-12-1.5 | Comb 12 Sample, 1.5mm |
| EL1400-12-2 | Comb 12 Sample, 2mm |

| Product Code | Description |
|------------------|---------------------------------|
| EL1400-14MC-0-75 | Comb & Load Guide 14/MC, 0.75mm |
| EL1400-14MC-1 | Comb 14 Sample MC, 1 mm |
| EL1400-14MC-1.5 | Comb & Load Guide 14/MC, 1.5mm |
| EL1400-14MC-2 | Comb 14 Sample MC 2mm |
| EL1400-16MC-0.75 | Comb 16 Sample MC, 0.75mm |
| EL1400-16MC-1 | Comb 16 Sample MC, 1mm |
| EL1400-16MC-1.5 | Comb 16 Sample MC, 1.5mm |
| EL1400-16MC-2 | Comb 16 Sample MC, 2 mm |
| EL1400-20-0.75 | Comb 20 Sample, 0.75mm |
| EL1400-20-1 | Comb 20 well, 1mm |
| EL1400-20-1.5 | Comb 20 Sample, 1.5mm |
| EL1400-20-2 | Comb 20 Sample, 2mm |
| EL1400-28MC-0-75 | Comb 28 well MC, 0.75mm |
| EL1400-28MC-1 | Comb 28 well MC, 1mm |
| EL1400-28MC-1.5 | Comb 28 Sample MC, 1.5mm |
| EL1400-28MC-2 | Comb 28 Sample MC, 2mm |
| EL1400-30MC-0.75 | Comb 30 Sample MC, 0.75mm |
| EL1400-30MC-1 | Comb 30 Sample MC, 1mm |
| EL1400-30MC-1.5 | Comb 30 Sample MC 1.5mm |
| EL1400-30MC-2 | Comb 30 Sample MC, 2mm |
| EL1400-35-0.75 | Comb 35 Sample MC, 0.75mm |
| EL1400-35-1 | Comb 35 Sample, 1mm |
| EL1400-35-1.5 | Comb 35 well, 1.5mm |
| EL1400-35-2 | Comb 35 Sample, 2mm |

Visit www.alphalabs.co.uk/electrophoresis to find out more.