

# TOP TIPS

for Cost -

## Four Practical Ways to Keep UN3373 Compliance Costs Down

- 1** Always use the smallest package possible to minimise mailing costs  
A one size fits all approach may seem like a good idea but you can be shipping a lot of air and inflating costs.



- 2** Obtain 95kPa certificates for primary containers  
Request these from your supplier to validate that your primary containers are 95kPa pressure resistance compliant. You can then select more cost-effective secondary containers.



- 3** Consider rigid secondary packaging  
The outer packaging then only needs to be an envelope. Lightweight polythene envelopes weigh a lot less than boxes and used with SpeciSafe® can minimise mailing costs.



- 4** Work with your transport provider  
Understand shipping costs before finalising packaging. Freight is often the largest portion of spend when transporting samples. To understand the total cost, freight charges and packaging costs need to be combined.



# Effective Sample Transport Packaging



## SpeciSafe® Secondary Mailing Packs

- Built-in absorbent reduces the number of components
- Only needs an outer envelope when used with 95kPa pressure tested vials
- Flat design minimises transport costs
- Fast and easy to use

