


MINI

VERTICAL GEL SYSTEMS




The Clarit-E Mini unit is perfect for routine protein electrophoresis. Sturdy construction, with simple, sliding gate clamps for gel casting. Runs up to 4 hand-cast or 2 pre-cast gels. Plates available with bonded spacers and colour coded with combs for compatibility.

KEY FEATURES:

The Clarit-E Mini provides a high-quality mini-gel electrophoresis system with a selection of interchangeable modules available for maximum versatility, enabling the tank to additionally run tube gels and perform electroblotting.

- α Available with 20 x 25cm, 20 x 20cm or 20 x 10cm gel trays
- α Run up to 550 samples
- α Low buffer volumes
- α Ideal for extended separations



Spacers and Combs

Colour coded by thickness, for easy identification and matching with combs

Systems and Modules:

Product Code	Description
EL2004	Gel Tank 10x10cm
EL2000	Dual Gel Tank 10x10cm + Casting Base
EL2012	Gel Tank 10x10cm, Casting Base + Ext Caster
EL2190	Quad System + Ext Caster
EL2400	Gel Tank 10x10cm System + Blotting Module
EL2320	Gel Tank 10x10cm System + 2D Tube Module

Product Code	Description
EL2200	Gel Tank System + 2D Tube & Blotting Modules
EL2404	Quad System + Blotting Module
EL2034	Blot Module Insert
EL2034-HI	Hi-Intensity Blot Module Insert
EL2028	2D Capillary IEF Gel Module Insert
EL2026	Gel Tank + 2D Capillary IEF Gel Module

Accessories:

Product Code	Description
EL2040	Mini Cooling Pack
EL2018	10x10cm Dual Caster
EL2080	Lid, No Cables
EL2082	Replacement Tank
EL2020	Replacement Caster Sealing Mats
EL2402	Blot Unit 10x10cm & 4 Cassettes
EL2036	Blotting Cassette
EL2038	Fibre Pads

Compatible Combs:

Product Code	Description
EL2000-10-1	Comb 10 Sample, 1mm
EL2000-10-1-5	Comb 10 Sample, 1.5mm
EL2000-1-1	Comb 1 Prep, 1mm
EL2000-12-0-5	Comb 12 Sample, 0.5mm
EL2000-12-1	Comb 12 Sample, 1mm
EL2000-16MC-1	Comb 16 well MC, 1mm
EL2000-16MC-1-5	Comb 16 well MC, 1.5mm
EL2000-20-0-5	Comb 20 Sample, 0.5mm
EL2000-20-1	Comb 20 well, 1mm
EL2000-10-1	Comb 10 Sample, 1mm
EL2000-12-1	Comb 12 Sample, 1mm
EL2000-20-1	Comb 20 well, 1mm
EL2000-1-1.5	Comb 1 Prep, 1.5mm
EL2000-5-1.5	Comb 5 Sample, 1.5mm
EL2000-8MC-1.5	Comb 8 well MC, 1.5mm
EL2000-9-1.5	Comb 9 Sample, 1.5mm
EL2000-10-1-5	Comb 10 Sample, 1.5mm
EL2000-12-1.5	Comb 12 Sample, 1.5mm
EL2000-20-1.5	Comb 20 Sample, 1.5mm
EL2000-1-2	Comb 1 Prep, 2mm
EL2000-5-2	Comb 5 Prep, 2mm
EL2000-8MC-2	Comb 8 well MC, 2mm
EL2000-9-2	Comb 9 Sample, 2mm
EL2000-10-2	Comb 10 Sample, 2mm
EL2000-12-2	Comb 12 Sample, 2mm
EL2000-20-2	Comb 20 Sample, 2mm

Clarit-E[®]

System and Power Supply Packages:

Product Code	Description
EL2078	Mini Vertical System with Power Supply
EL2430	Mini Vertical System + Blotting Module with Power Supply

Compatible Power Supplies and Adapters:

Product Code	Description
EL2707	Mini Power Supply 500V, 400mA, 120W
EL2705	Mini-300 Power Supply 300V, 400mA, 60W
EL2702	Power Pack Adapters 4mm to 2mm
EL2752	Maxi Power Pack 300V, 3A, 300W
EL2732	Midi Power Pack 300V, 700mA, 150W
EL2714	Replacement Cables

Glass Plates and Gel Tubes:

Product Code	Description
EL2042	10x10cm Notched Glass Plates 2mm
EL2044	10x10cm Plain Glass Plates 2mm
EL2046	10x10cm Plain Plates + 0.5mm Spacers
EL2048	10x10cm Plain Plates + 0.75mm Spacers
EL2050	10x10cm Plain Plates + 1mm Spacers
EL2052	10x10cm Plain Plates + 1.5mm Spacers
EL2054	10x10cm Plain Plates + 2mm Spacers
EL2064	10x10cm Notched Plates + 2mm Spacers
EL2062	10x10cm Notched Plates + 1.5mm Spacers
EL2058	10x10cm Notched Plates + 0.75mm Spacers
EL2056	10x10cm Notched Plates + 0.5mm Spacers
EL2060	10x10cm Notched Plates + 1mm Spacers
EL2066	Dummy Plate, 10x10cm
EL2068	0.5mm Spacers for 10cm Glass Plates
EL2070	0.75mm Spacers for 10cm Glass Plates
EL2072	1mm Spacers for 10cm Glass Plates
EL2074	1.5mm Spacers for 10cm Glass Plates
EL2076	2mm Spacers for 10cm Glass Plates
EL2030	Mini Capillary Tubes 1mm Internal Diam
EL2356	Mini Capillary Tubes 1.5mm Internal Diam
EL2032	Capillary Blanking Ports

Visit www.alphalabs.co.uk/electrophoresis to find out more.