


## CHOICE

## HORIZONTAL GEL SYSTEMS



The Clarit-E Choice offers a wide degree of versatility, with three tray options available allowing the choice of one, two or all three gel length options. The 15cm total run length allows restriction fragment or other close MW sample bands to be easily separated and identified.

### KEY FEATURES:

The Clarit-E Choice is ideal for restriction fragment analysis, sample prep or checking of high numbers of samples.

- α Three tray options
- α Run up to 210 samples
- α Low buffer volumes
- α Multichannel pipette compatible combs for speed loading



**Casting Dams**  
Cast gels separately from the gel tank and save time while the tank is in use.



**Buffer Saver Blocks**  
Save costs by physically reducing the volume of buffer required in the gel chamber.

### Tanks:

Product Code	Description
EL1400	Gel Tank 15x10cm UV Tray
EL1450	Gel Tank 15x7cm UV Tray
EL1500	Gel Tank 15x15cm UV Tray
EL1510	Gel Tank Trio, with 3x 15cm Tray

### Trays:

Product Code	Description
EL1402	UV Tray 15x10cm
EL1452	UV Tray 15x7cm
EL1502	UV Tray 15x15cm


## Accessories:

Product Code	Description
EL1404	Casting Dams
EL1406	Electrode, positive
EL1408	Electrode, negative
EL1410	Loading Guides
EL1412	Viewing Platform
EL1414	Cool Pack and Platform
EL1416	Buffer Saver Blocks
EL1420	UV Gel Scoop
EL1418	Flexi Caster
EL1026	Gel Levelling Table
EL1422	Lid

# Clarit-E®

## Compatible Power Supplies and Adapters:

Product Code	Description
EL2707	Mini Power Supply 500V, 400mA, 120W
EL2705	Mini-300 Power Supply, 300V, 400mA, 60W
EL2702	Power Pack Adapters 4mm To 2mm
EL2752	Maxi Power Pack 300V, 3A, 300W
EL2732	Midi Power Pack 300V, 700mA, 150W
EL2714	Replacement Cables



## Compatible Combs:

Product Code	Description
EL1400-1-0.75	Comb Prep 1, Marker 1, 0.75 mm
EL1400-1-1	Comb Prep 1, Marker 1, 1 mm
EL1400-1-1.5	Comb Prep 1, Marker 1, 1.5mm
EL1400-1-2	Comb Prep 1, Marker 1, 2mm
EL1400-2-0.75	Comb Prep 2, Marker 2, 0.75mm
EL1400-2-1	Comb Prep 2, 1.0mm
EL1400-2-1.5	Comb Prep 2, Marker 2, 1.5mm
EL1400-2-2	Comb Prep 2, Marker 2, 2mm
EL1400-4-0.75	Comb Prep 4, Marker 2, 0.75 mm
EL1400-4-1	Comb Prep 4, Marker 2, 1 mm
EL1400-4-1-5	Comb Prep 4, Marker 2, 1.5 mm
EL1400-4-2	Comb Prep 4, Marker 2, 2 mm
EL1400-10-0.75	Comb 10 Sample, 0.75mm
EL1400-10-1	Comb 10 well, 1 mm
EL1400-10-1.5	Comb 10 well, 1.5 mm

Product Code	Description
EL1400-10-2	Comb 10 well, 2 mm
EL1400-10MC-0.75	Comb 10 Sample MC, 2 mm
EL1400-10MC-1	Comb 10 well MC, 1mm
EL1400-10MC-1.5	Comb 10 well MC, 1.5mm
EL1400-10MC-2	Comb 12 well MC, 0.75mm
EL1400-12-0.75	Comb 10 well MC, 1.5mm
EL1400-12-1	Comb 12 Sample, 1mm
EL1400-12-1.5	Comb 12 Sample, 1.5mm
EL1400-12-2	Comb 12 Sample, 2mm
EL1400-14MC-0-75	Comb & Load Guide 14/MC, 0.75mm
EL1400-14MC-1	Comb 14 Sample MC, 1 mm
EL1400-14MC-1.5	Comb & Load Guide 14/MC, 1.5mm
EL1400-14MC-2	Comb 14 Sample MC, 2mm
EL1400-16MC-0.75	Comb 16 Sample MC, 0.75mm
EL1400-16MC-1	Comb 16 Sample MC, 1mm

Visit [www.alphalabs.co.uk/electrophoresis](http://www.alphalabs.co.uk/electrophoresis) to find out more.