


MINI WIDE VERTICAL GEL SYSTEMS



With a generous 20cm gel width, the Clarit-E Mini Wide vertical gel unit enables resolution of up to 48 samples per gel, ideal for higher throughput applications. For maximum functionality, interchangeable inserts are offered to enable electroblotting or 2D gel tube electrophoresis using the same base tank.

KEY FEATURES:

The Clarit-E Mini Wide uses 4mm thick glass plates to reduce breakage with bonded spacers for ease of assembly and casting.

- α The capacity of two standard mini gels in a single gel
- α Compact system with economical buffer usage
- α Interchangeable insert modules for versatility



Insert Modules

One tank, three applications – PAGE, 2D
Tube Gels, Electroblotting

Systems and Modules:

Product Code	Description
EL2083	Gel Tank 20x10cm
EL2084	Mini Wide System 20x10cm + Casting Base
EL2086	Mini Wide System 20x10cm, Casting Base + Ext Caster
EL2450	Mini Wide System 20x10cm + Blotting Module
EL2330	Mini Wide System 20x10cm + 2D Tube Module
EL2250	Mini Wide System + 2D Tube & Blotting Modules
EL2114	Blot Module Insert
EL2098	2D Capillary IEF Gel Module Insert
EL2096	Gel Tank with 2D Capillary IEF Gel Module


Glass Plates:

Product Code	Description
EL2112	20x10cm Notched Glass Plates
EL2116	20x10cm Plain Glass Plates
EL2118	20x10cm Plain Glass Plates + 0.75mm Spacers
EL2119	20x10cm Plain Glass Plates + 1mm Spacers
EL2122	20x10cm Plain Plates + 1.5mm Spacers
EL2125	20x10cm Plain Plates + 2mm Spacers
EL2127	20x10cm Notched Plates + 0.75mm Spacers
EL2129	20x10cm Notched Plates + 1mm Spacers
EL2156	20x10cm Dummy Plate 4mm

Accessories:

Product Code	Description
EL2068	0.5mm Spacers for 10cm Glass Plates
EL2070	0.75mm Spacers for 10cm Glass Plates
EL2072	1mm Spacers for 10cm Glass Plates
EL2074	1.5mm Spacers for 10cm Glass Plates
EL2076	2mm Spacers for 10cm Glass Plates
EL2018	10cm Dual Caster
EL2094	Inner Run Module
EL2106	Blotting Cassette
EL2030	Mini Capillary Tubes 1mm Internal Diam
EL2356	Mini Capillary Tubes 1.5mm Internal Diam
EL2032	Capillary Blanking Ports
EL2452	Blot Unit 20x10cm & 4 Cassettes
EL2106	Blotting Cassette
EL2458	Fibre Pads

Clarit-E®



Compatible Combs:

Product Code	Description
EL2100-18MC-1	Comb 18 well MC, 1mm
EL2100-18MC-1.5	Comb 18 well MC, 1.5mm
EL2100-18MC-2	Comb 18 well MC, 2mm
EL2100-24-1	Comb 24 Sample, 1mm
EL2100-24-1.5	Comb 24 Sample, 1.5mm
EL2100-30-1	Comb 30 Sample, 1mm
EL2100-36-1-5	Comb 36 Sample, 1.5mm
EL2100-36MC-1	Comb 18 well MC, 1mm
EL2100-48-1.5	Comb 36 Sample, 1.5mm

Compatible Power Supplies and Adapters:

Product Code	Description
EL2707	Mini Power Supply 500V, 400mA, 120W
EL2705	Mini-300 Power Supply 300V, 400mA, 60W
EL2702	Power Pack Adapters 4mm to 2mm
EL2752	Maxi Power Pack 300V, 3A, 300W
EL2732	Midi Power Pack 300V, 700mA, 150W
EL2714	Replacement Cables

Visit www.alphalabs.co.uk/electrophoresis to find out more.