

PRODUCT DATA SHEET

GRADUATED PASTETTE®

GRADUATED TRANSFER PIPETTES IN A RANGE OF SIZES FOR QUICK AND EASY MEASURED TRANSFERS OF FLUIDS AND REAGENTS.



Ideal for:

- Cytology and slide staining techniques
- Measured transfers in microbiology
- Clinical sampling from blood bag tubing
- Sample preparation in chromatography

Graduated Pastettes® are the easiest and most cost effective way of making quick and convenient measured transfers.

High quality manufacturing from virgin low-density polyethylene means graduated Pastettes are versatile and safe for you to use with a wide variety of liquids.

Fractional delivery of small volumes is achieved with the soft responsive bulb allowing drop-wise control of measured volumes.



Frequently referred to in National Standard Methods and NHS SOP, sterile and non-sterile Pastettes are used in many laboratory and diagnostic procedures as an effective solution for specific liquid handling tasks. Graduated Pastettes are available in a choice of forms in volumes of 1, 2, 3, & 5mls. See over for the complete range.

Sterile Packs

Sterile packs are securely sealed and e-beam irradiated to provide a sterility assurance level of 10⁻⁴.


Pack formats of single wrapped, 10's and 20's, provide the ideal solution to aseptic liquid handling procedures where sample protection is of key importance.

Certificates of sterility are available by lot number.

A choice of sterile pack sizes provides cost effective quantities to suit your aseptic procedures.

See reverse for full range of Graduated Pastettes

	1ml Graduated	3ml Graduated	5ml Graduated	Graduated Medi	Small Bulb Graduated	Graduated Mini
Length (mm)	155	155	217	152	140	116
Drops per ml	25	20	19	24	25	21
Graduations at	0.25ml	0.50ml	1.00ml	0.50ml	0.25ml	0.1ml



Are you using genuine Pastettes®?

Often imitated, but never bettered, Alpha Laboratories offer the most comprehensive range of disposable transfer pipettes available to suit all clinical, experimental and routine liquid handling procedures.

Pastettes are manufactured from virgin low-density polyethylene (LDPE), safely providing precise transfer and delivery of liquids, reagents and cell suspensions. Graduated, Non-graduated, Fine Tip and Micro ranges feature designs specific for your applications.

Pastettes are ideally suited to sterile and non-sterile applications. Available in sterile packs of 10's, 20's or single wrapped as well as non-sterile bulk packaging.

Request the Alpha Laboratories' catalogue to see our complete range or request samples to try before you buy.

Non Sterile						
Cat No	LW4040	LW4111	LW4728	LW4274	LW4212	LW4282
Pack Size / IRLA Price	500	500	250	500	500	500
Sterile 1's						
Cat No	LW4005	LW4112				LW4284
Pack Size / IRLA Price	500	500				500
Sterile 10's						
Cat No	LW4006	LW4113				
Pack Size / IRLA Price	500	500				
Sterile 20's						
Cat No	LW4010	LW4114	LW4729	LW4277	LW4215	LW4283
Pack Size / IRLA Price	500	500	240	500	500	500

Non-sterile

1's

10's

20's