

# CD-Chex™ Daily

One vial. 51 markers. Your panel covered.



## Control Marker Guide/Coverage

| Phenotype         | CD-Chex® Daily* | CD-Chex Plus® | CD-Chex CD34® | CD-Chex CD103 Plus** | CD-Chex CD117® Plus* | CD-Chex TdT Plus® |
|-------------------|-----------------|---------------|---------------|----------------------|----------------------|-------------------|
| CD1a              |                 |               |               |                      |                      | •                 |
| CD2               | •               | •             |               |                      |                      | •                 |
| CD3               | •               | •             |               |                      |                      |                   |
| cCD3              |                 |               |               |                      |                      | •                 |
| CD4               | •               | •             |               |                      |                      | •                 |
| CD5               | •               | •             |               |                      |                      | •                 |
| CD7               | •               | •             |               |                      |                      | •                 |
| CD8               | •               | •             |               |                      |                      | •                 |
| CD10              | •               | •             |               |                      |                      |                   |
| CD11b             | •               | •             |               |                      |                      |                   |
| CD11c             | •               | •             |               |                      |                      |                   |
| CD13              | •               | •             |               |                      |                      |                   |
| CD14              | •               | •             |               |                      |                      |                   |
| CD15              | •               | •             |               |                      |                      |                   |
| CD16              | •               | •             |               |                      |                      |                   |
| CD16-CD56         | •               | •             |               |                      |                      |                   |
| CD19              | •               | •             |               |                      |                      |                   |
| CD20              | •               | •             |               |                      |                      |                   |
| CD22              | •               | •             |               |                      |                      |                   |
| CD23              | •               | •             |               |                      |                      |                   |
| CD24              | •               |               |               |                      |                      |                   |
| CD25              |                 |               |               |                      | •                    |                   |
| CD26              | •               |               |               |                      |                      |                   |
| CD27              | •               |               |               |                      |                      |                   |
| CD28              | •               |               |               |                      |                      |                   |
| CD30              |                 |               |               | •                    |                      |                   |
| CD33              | •               | •             |               |                      |                      |                   |
| CD34 (Progenitor) | •               | •             | •             |                      |                      |                   |
| CD34 (Abnormal)   |                 |               |               |                      |                      | •                 |
| CD36              | •               |               |               |                      |                      |                   |
| CD38+ / CD45+     | •               | •             |               | •                    |                      |                   |
| CD38+ / CD45-     |                 |               |               | •                    |                      |                   |
| CD41              | •               |               |               |                      |                      |                   |
| CD43              | •               |               |               |                      |                      |                   |
| CD45RA            | •               |               |               |                      |                      |                   |
| CD45RO            | •               |               |               |                      |                      |                   |
| CD45              | •               | •             |               |                      |                      |                   |
| CD56              | •               | •             |               | •                    |                      |                   |
| CD57              | •               |               |               |                      |                      |                   |
| CD58              | •               |               |               |                      |                      |                   |
| CD61              | •               |               |               |                      |                      |                   |
| CD64              | •               |               |               |                      |                      |                   |
| CD71              |                 |               |               |                      | •                    |                   |
| CD79a             | •               |               |               |                      |                      |                   |
| CD79b             | •               |               |               |                      |                      |                   |
| CD81              | •               |               |               |                      |                      |                   |
| CD103             |                 |               |               | •                    |                      |                   |
| CD117             |                 |               |               |                      | •                    |                   |
| CD138             |                 |               |               | •                    |                      |                   |
| CD200             | •               |               |               |                      |                      |                   |
| FMC7              | •               |               |               |                      |                      |                   |
| HLA-DR            | •               | •             |               |                      |                      |                   |
| IgM               | •               |               |               |                      |                      |                   |
| Kappa / CD19      | •               | •             |               |                      |                      |                   |
| Lambda / CD19     | •               | •             |               |                      |                      |                   |
| cKappa            | •               |               |               |                      |                      |                   |
| cLambda           | •               |               |               | •                    |                      |                   |
| cMPO              | •               |               |               |                      |                      |                   |
| TCR α/β           | •               |               |               |                      |                      |                   |
| TCR γ/δ           | •               |               |               |                      |                      |                   |
| TdT               |                 |               |               |                      |                      | •                 |
| TRBC1             | •               |               |               |                      |                      |                   |
| Lymphocyte        | •               | •             |               |                      |                      |                   |
| Monocyte          | •               | •             |               |                      |                      |                   |
| Granulocyte       | •               | •             |               |                      |                      |                   |
| WBC               | •               | •             | •             | •                    | •                    | •                 |

• CD-Chex Plus RUO markers will be discontinued with lot 6348, shipping Dec. 14, 2026.

\*CD-Chex Daily, CD-Chex Select, CD-Chex CD103 Plus, CD-Chex CD117 Plus, and CD-Chex TdT Plus are For Research Use Only. Not for use in diagnostic procedures.