



Specimen Containers

For the safe collection, storage and transportation of biological samples

- α Specially developed for biological sample collection and processing
- α Ideal for applications where strict adherence to regulatory requirements is essential
- α Leak-proof tested to IATA 95kPa
- α Food contact approved
- α CE IVD marked



Widely used in a range of applications including:

- Pathology
- Histology
- Water Sampling
- Metals Analysis

Sturdy Design

The Specimen Containers are made from durable polypropylene with graduations in 10ml increments on 3 sides, to ensure contents and their volumes are easily visible.

The caps are manufactured from high density polyethylene (HDPE).

Both container and cap have easy-to-grip ridges to ensure opening and closing the container is easy, even when wearing gloves.

These robust containers are ideal for liquid or solid biological specimen collection, storage and transportation in harsh conditions.

Testing and Validation

The Specimen Containers are food contact approved and testing has been carried out for suitability with aggressive chemicals such as formaldehyde.

A process validation has been carried out on these disposable plastic specimen containers to ensure the unique "double seal" closure technology remains leak-proof under a pressure differential of not less than 95 kPa (0.95 bar, 14 psi).

The containers together with their caps are certified to be leak-proof to IATA 95kPa.

This certification makes them suitable for sample transport according to UN3373 packaging instructions, when combined with appropriate secondary and tertiary packaging, e.g. SpeciSafe® or other Rigid Flexible Packaging Systems.

Specifications

COMPATIBILITY: Specimen Container Caps are only compatible with these Specimen Containers.

VOLUME: 120ml

CONTAINER HEIGHT WITHOUT CAP: 75mm

HEIGHT WITH CAP: 77mm

DIAMETER WITHOUT CAP: 54mm

DIAMETER WITH CAP: 57mm

MATERIAL: Polypropylene container and HDPE cap

COLOUR: Semi-transparent Container
White Cap

STERILITY: Non Sterile

LEAK-PROOF: Tested to IATA 95kPa

CE IVD MARK: Yes

PACK SIZE: Supplied as 500 containers or caps per pack.
Containers and caps sold separately.

Ordering Information

| Cat. No. | Description | Pack Size |
|-----------|--|-----------|
| STC0120CW | Specimen Container Cap 120ml White Non-Sterile | 500 |
| STC0120SJ | Specimen Container Jar 120ml Clear Non-Sterile | 500 |



Find out more at www.alphalabs.co.uk

Manufactured for Alpha Laboratories by Splice Cast Ltd, Unit F, Mochdre Industrial Estate, Newtown Powys, SY16 4LE.