



## Custom Designed Consumables

*Specialist Components for Clinical Trials,  
Diagnostic and Laboratory Kit Manufacture*

You define:

- Specifications
- Size
- Branding
- Packaging

*Find out how we can help you with a range of bespoke products*



## Disposable Plastic Liquid Transfer Pipettes

The Alpha brand Pastette® is the leading name for disposable liquid transfer or Pasteur pipettes. It is available in over 50 standard designs, many developed in response to specific customer requirements.

The Dual-Bulb Pastette range is just one example. It is the ideal pipette for inclusion in assay kits or for use in other applications requiring fast and easy delivery of an exact measured volume. This simple solution can readily be used in untrained hands and is available in a range of volumes from 20 to 300µl.

In addition Alpha Laboratories offers a customised design service to develop pipettes specifically to OEM requirements.



## UN3373 Compliant Sample Transport Packaging

Our SpecSafe® mailing packs, and other packaging solutions can be developed into custom designed kits that are ideal for clinical trials and diagnostic assays. They ensure convenience and complete compliance for transporting biological samples between participants and testing laboratories.



## Plate Sealing Solutions

Diagnostic assay kits often contain reagents held in microtitre or PCR plates that need to be securely sealed.

Alpha Laboratories provides a range of plate sealing options that can be customised for specific manufacturing and packaging requirements. These include indexed roll formats for automation, speciality materials and adhesives, custom die-cut shapes and sizes, plus custom printing, packaging and labelling.

Stringent QC procedures for all products and certificates of compliance and sterility provide assurance of product consistency lot-to-lot.



## Reagent Tubes, PCR and Other Consumables

We also offer a wide range of screw cap microtubes, PCR plates and tubes, Universal and Bijou containers etc. that may also be required for kit manufacture.

## Customised Sample Transport Packs for Clinical Trials and Contract Research Projects

The regulatory requirements for shipping clinical trial samples can sometimes be troublesome and a time-consuming distraction from the key activity of developing life-saving drugs.

However, whether you are involved in National or Global projects, transporting clinical trial samples safely, securely and in adherence to current regulations is essential.

Alpha Laboratories can ease this burden and deliver solutions which are unique, compliant with transport regulations, easy to use and bespoke to your requirements – allowing you to spend more time on your most important tasks.

## Bespoke UN3373 Compliant Transport Packaging for Biological Specimens

The innovative SpecSafe® all in one secondary packaging system combines absorbent material, vial separation and rigid container within a single pre-moulded mailing pack. SpecSafe packs are available in a range of designs to fit perfectly around commonly used 95 kPa sample vessels. They must be used with primary containers which are 95kPa pressure differential certified. A rigid outer container is required to complete the UN3373 compliant solution.

SpeciSafe® packs can be used by Clinical Trial Central Laboratories, Diagnostics Reference Laboratories, CROs, plus Pharmaceutical and Biotech companies for their clinical kit preparation and regular distribution needs.

They can also be supplied in custom designed packs as a supply chain logistics solution for diagnostics, clinical trials and contract research projects in the pharmaceutical, healthcare and laboratory industries.





## UN3373 Compliant Secondary Packaging

Alpha Laboratories also offers a wide range of other UN3373 compliant secondary packaging such as rigid boxes, absorbent sheets, 95kPa pouches, labels, mailing bags and envelopes.

These custom designed transport packaging materials are also ideal for use as components in diagnostic or laboratory kit manufacture.

- SpeciSafe® UN3373 Compliant Mailing Packs
- Complete Mailing Systems
- Absorbent Sheets
- 95kPa Pouches
- Rigid Outer Boxes
- Labels and Seals
- Mailing Bags and Envelopes
- Medical Specimen Collection Containers
- Biohazard Bags
- Specimen Bags
- Slide Mailers
- Transport Tubes

**Contact Alpha Laboratories for a quotation for your custom requirements and join our long list of satisfied customers.**

### Alpha Laboratories Limited

40 Parham Drive  
Eastleigh, Hampshire  
SO50 4NU United Kingdom

**Tel:** +44 (0) 23 8048 3000

**Fax:** +44 (0) 23 8064 3701

**Email:** [sales@alphalabs.co.uk](mailto:sales@alphalabs.co.uk)

**Web:** [www.alphalabs.co.uk](http://www.alphalabs.co.uk)

### Customer Account Number

|  |  |  |  |  |  |
|--|--|--|--|--|--|
|  |  |  |  |  |  |
|--|--|--|--|--|--|