



**TECHNOPATH**  
CLINICAL DIAGNOSTICS

**Multichem<sup>®</sup> ID**  
INFECTIOUS DISEASE QC

POWERED BY



# Multichem<sup>®</sup> ID-SeroNeg

Providing Instrument Specific Third-Party  
Infectious Disease Quality Controls



Available in Europe.

## Features

- ◆ Ready to use liquid stable controls
- ◆ Direct sampling from control vials
- ◆ Instrument specific barcodes
- ◆ Technical and scientific support

## Specifications

- ◆ 36 month closed vial stability at 2°C to 8°C.
- ◆ 30 day open vial stability at 2°C to 8°C.

Product	Description	Config.	Part Code
Multichem ID-SeroNeg	Negative Control	4 x 5mL	SR100N
Multichem ID-SeroNeg Mini Kit	Negative Control Mini Kit	1 x 5mL	SR100MN

For further information see Instructions For Use.

*THE* QUALITY CONTROL COMPANY

# Multichem<sup>®</sup> ID-SeroNeg

## Analyte List

### Analytes

---

Anti-HIV

Anti-HCV

Anti-HBc

Anti-HTLV

Anti-Treponema

HBsAg

HIV p24



**TECHNOPATH**  
CLINICAL DIAGNOSTICS

TECHNOPATH CLINICAL DIAGNOSTICS  
Technopath Life Sciences Park, Fort Henry, Ballina, Co. Tipperary, V94 FF1P, Ireland.  
For Multichem QC Support: [qcsupport@technopathcd.com](mailto:qcsupport@technopathcd.com) Website: [www.technopathcd.com](http://www.technopathcd.com)