


The next generation in QC Data Management from Technopath Clinical Diagnostics

IAMQC® Infinity is a powerful **secure** cloud based QC Data management solution.

Our 'software as a service' (SAAS) model ensures an experience that is:

- Light on upfront capital expenditure in the lab
- Light on lab IT resource
- Instant access to latest releases and updates



Consolidating all QC data management functions into one scalable solution

IAMQC Infinity is designed as an open QC data management solution that complements and supports Technopath Clinical Diagnostics Multichem® Quality Control (QC) product range, and provides Laboratory Managers and Technologists with a range of QC software tools to analyse their QC results in real-time.



Feature & Functions


- Peer comparison
- Daily QC management
- Westgard rules
- QC reports
- Six Sigma analysis
- Total Error
- Flexible reports
- QC Scheduler
- Recommended rules wizard
- Reagent verification

IAMQC Infinity is practical, graphical, user-definable and easy to use and provides laboratories with a centralised standardised approach to QC data management.

Your QC Control Central!

The most important requirement in managing QC data is knowing immediately when and where you have an issue. IAMQC Infinity brings this to a new level with a dashboard showing a traffic light system of performance for each instrument and highlighting what tests have issues. There is also a notification that will appear showing a consolidated list of all instruments and tests with issues, allowing those issues to be prioritised and dealt with in a timely fashion.

- Dashboard view to show overall performance across multiple labs – highlighting issues quickly and clearly
- Flexible reports engine that allows users to pick the data relevant to them and how the data should be displayed – by instrument, by test, by lot number etc.
- Easy access reports to desired output – either PDF or Microsoft Excel




IAMQC Infinity gives the end user multiple layers of access to information by providing a high level overview of performance at the dashboard level and guiding the user into more details at the instrument and test level.

IAMQC infinity set up could not be simpler!

IAMQC Infinity uses a simple registration form capturing minimal details from the customer in order to get activated on the system. Once the information is submitted and the registration has been activated – you're ready to use the system!

A wizard based setup form guides the lab user through entering all the required information to get your lab configured. Simple tabs along the top of the wizard show you the area you are configuring and allows you to move seamlessly through each area.

- Easy to use intuitive web based system
- Administration area to hide access to functions not required by specific users.
- Easy to deploy, simple registration form.
- Straightforward setup wizard to guide you through the system requirements.
- Import setup template wizard to allow for importation of instruments and products.



IAMQC Infinity is a trade mark of Technopath Clinical Diagnostics.