



TECHNOPATH
CLINICAL DIAGNOSTICS

Multichem® ID
INFECTIOUS DISEASE QC

POWERED BY



Multichem® ID-B

Providing Instrument Specific Third-Party
Infectious Disease Quality Controls



Available in Europe.

Features

- ◆ Ready to use liquid stable controls
- ◆ Direct sampling from control vials
- ◆ Instrument specific barcodes
- ◆ Technical and scientific support

Specifications

- ◆ 24 month closed vial stability at 2°C to 8°C.
- ◆ 30 day open vial stability at 2°C to 8°C.

Product	Description	Config.	Part Code
Multichem ID-B	Positive Control	4 x 5mL	SR101B
Multichem ID-B Mini Kit	Positive Control Mini Kit	1 x 5mL	SR101MB

For further information see Instructions For Use.

THE QUALITY CONTROL COMPANY

Multichem[®] ID-B

Analyte List

Analytes

Anti-HIV 1
 Anti-HCV
 Anti-HBc
 HBsAg
 Anti-HTLV
 Anti-Treponema Pallidum

Analyser Systems

MANUFACTURER	SYSTEM	ASSAY NAME
Abbott	ALINITY i	Anti-HBc II
		Anti-HCV
		HBsAg Qualitative II
		HIV Ag/Ab Combo
		rHTLV I/II
		Syphilis TP
Abbott	ALINITY s	Anti-HBc
		Anti-HCV
		HBsAg
		HIV Ag/Ab Combo
		HTLV I/II
		Syphilis
Abbott	ARCHITECT	Anti-HBc II
		Anti-HCV
		HBsAg Qualitative II
		HIV Ag/Ab Combo
		rHTLV-I/II
		Syphilis TP
DiaSorin	LIAISON	Anti-HBc
		Treponema Screen
DiaSorin	LIAISON XL MUREX	HBsAg Quant
		HCV Ab
		HIV Ab/Ag
		recHTLV-I/II



TECHNOPATH
 CLINICAL DIAGNOSTICS

TECHNOPATH CLINICAL DIAGNOSTICS

Technopath Life Sciences Park, Fort Henry, Ballina, Co. Tipperary, V94 FF1P, Ireland.
 For Multichem QC Support: qcsupport@technopathcd.com Website: www.technopathcd.com